



SPECIFICATION SHEET Foam Cleaner

Professional Lysol® Foam Cleaner - Aerosol

Foam cleaning action you can see working. Clings to vertical surfaces.



Product Information

Product	Units/Case	Unit Size	Reckitt Item #
Professional Lysol® Foam Cleaner - Aerosol	12	24 oz	36241-02775

Shipping Specifications

Reckitt Item #	36241-02775
Package Size	24 oz
Case Pack	12
Pallet Cases	85
Tie/High	17/5
Case Weight	23.40 lbs
Case Cube	0.64 cu. ft
Case L x W x H	11.88 x 8.75 x 10.63 in
Item D x W x H	2.81 x 2.81 x 10.13 in
Item Cube	0.046 cu. ft
Item Weight	1.89 lbs
Item UPC	0-36241-02775-9
CSC Case GTIN	10036241027756

Technical Specifications

Technical Data	Professional Lysol® Foam Cleaner - Aerosol
Color/Form	Clear/Liquid
Fragrance	Fragrant
Active Ingredient	Didecyl dimethyl ammonium chloride (0.025%) Didoctyl dimethyl ammonium chloride (0.025%) Octyl decyl dimethyl ammonium chloride (0.050%)
Other Ingredients	99.900%
pH	12.2 to 13.2
EPA Registration	777-71-675

DIRECTIONS FOR USE¹

Pre-clean. Apply and leave wet to sanitize² and 5 minutes to disinfect. Wipe dry. See list of organisms for appropriate contact times.

USE AREAS	SURFACE TYPES ³
<ul style="list-style-type: none"> - Appliance exteriors - Bathtubs & showers - Cabinets - Chairs - Bathroom counters - Door knobs - Exhaust fans - Floors - Garbage cans 	<ul style="list-style-type: none"> - Light fixtures - Stove tops - Tables - Toilet bowl exteriors - Urinals - Vinyl Shower Curtains - Windows - Other, hard, non-porous surfaces
	<ul style="list-style-type: none"> - Chrome Fixtures - Crystal - Enamel - Fiberglass fixtures - Glass - Glazed Ceramic/ Porcelain Tile - Hard Plastic - Linoleum - Metal - Synthetic Marble - Chrome/Stainless Steel Fixtures - Vinyl - Non-Medical Chrome/ Stainless Steel

¹On hard, non-porous surfaces

²See list of organisms for appropriate contact times.

GERMICIDAL EFFICACY

5 Minute Disinfection - Bacteria	Ref. No.
Pseudomonas aeruginosa	ATCC 15442
Salmonella enterica	ATCC 10708
10 Minute Disinfection - Bacteria	Ref. No.
Acinetobacter baumannii Multi-Drug Resistant (MDR)	ATCC 19606
Bordetella pertussis	ATCC 12743
Enterobacter aerogenes Multi-Drug Resistant (MDR)	ATCC 29751
Enterococcus faecalis	ATCC 828
Enterococcus faecium Multi Drug Resistant (MDR)	ATCC 51559
Escherichia coli O157:H7	ATCC 43888
Escherichia coli (ESBL) (Extended- Spectrum beta-lactamase) (producing) (Escherichia coli)	ATCC BAA-196
Klebsiella pneumoniae (New Delhi metallo- beta-lactamase 1) (NDM-1)	CDC 1000527
Listeria monocytogenes	TCC 7644
Staphylococcus aureus	ATCC 6538
(Methicillin Resistant) Staphylococcus aureus (MRSA)	ATCC 33592
(Community Acquired Methicillin Resistant) Staphylococcus aureus (CA-MRSA) (Genotype USA 300)	NARSA NRS 384
Vancomycin Resistant Staphylococcus aureus (VRSA)	NARSA VR51
Streptococcus pyogenes	ATCC 12384
10 Second Sanitization - Bacteria (Hard, Non-Porous, Non-Food Contact Surfaces)	Ref. No.
Enterobacter aerogenes	ATCC 13048
Escherichia coli	ATCC 11229
Staphylococcus aureus	ATCC 6538
30 Second Sanitization - Bacteria (Hard, Non-Porous, Non-Food Contact Surfaces)	Ref. No.
(Vancomycin Resistant Enterococcus) Enterococcus faecalis (VRE)	ATCC 51299
Escherichia coli O157:H7	ATCC 43888
(Methicillin Resistant) Staphylococcus aureus (MRSA)	ATCC 33592
10 Minute Disinfection - Fungi	Ref. No.
Candida albicans	ATCC 10231
Hard Surface Mildew Fungistat	Ref. No.
Aspergillus niger (mold & mildew)	ATCC 6275

5 Minute Disinfection - Viruses	Ref. No.
Bovine Viral Diarrhea Virus (Strain Oregon C24v-genotype 1) as a Surrogate Virus for Human Hepatitis C Virus	N/A
2009-H1N1 Influenza A Virus (Novel H1N1) (Strain A/Mexico/4108/2009)	CDC 2009712192
Feline Infectious Peritonitis Virus (Strain FIPV-T146, P100)	ATCC VR-2202
Herpes Simplex Virus Type 1 (Strain F (I))	ATCC VR-733
Herpes Simplex Virus Type 2 (Strain G)	ATCC VR-734
Influenza A Virus (Stain Hong Kong)	ATCC VR-544
Influenza A Virus (Strain A/New Caledonia/20/99)	CDC 99044140
Influenza B Virus (Strain B/Hong Kong/5/72)	ATCC VR-823
Avian Influenza A Virus (H5N1) (Strain VN5N1-PR8/CDC-RG)	CDC 2006719965
Avian Influenza A (H3N2) Reassortant Virus (Stain A/Washington/897/80 x A/Mallard/New York/6750/78)	ATCC VR-2072
Human Immunodeficiency Virus Type 1 (Strain HTLV-III B)	N/A
Respiratory Syncytial Virus (RSV) (Strain Long)	ATCC VR-26
(Severe Acute Respiratory Syndrome) SARS-Associated Coronavirus	CDC 200300592
10 Minute Disinfection - Viruses	Ref. No.
Adenovirus Type 2 (Strain Adenoid 6)	ATCC VR-846
Canine Parvovirus (Strain Cornell - 780916-80)	ATCC VR-2017
Coxsackievirus Type B3 (Strain Nancy)	ATCC VR-30
Echovirus Type 12 (Strain Travis 2-85)	ATCC VR-1563
Feline Calicivirus (as a surrogate virus for) (Norovirus)	ATCC VR-782
Feline Panleukopenia Virus (Strain Philips-Roxane)	ATCC VR-648
Feline Picornavirus (Strain FRV (17FRV))	ATCC VR-649
Hepatitis A Virus (Strain HM-175)	N/A
Poliovirus Type 1 (Strain Chat)	ATCC VR-1562
Porcine Parvovirus (PPV) (Strain NADL-2)	ATCC VR-742
Rhinovirus Type 39 (Strain 209)	ATCC VR-340
Rotavirus (Strain) WA	ATCC VR-2018
Duck Hepatitis B Virus as a Surrogate Virus for Human Hepatitis B Virus	N/A

For additional product information and Safety Data Sheet, please visit www.rbnainfo.com

